



Application of MALDI-FT-ICR-MS for nanostructured thin films characterization

EU_FT-ICR-MS 2nd Advanced Users School

Bogdan Purcareanu¹, Alina Morosan², Adrian Costache², Alexia Popescu², Bogdan Mihaiescu³, Rodica Cristescu⁴, Daniela Istrati², Dan Eduard Mihaiescu²

1 - Biotehnos SA, Otopeni city,

2 - University Politehnica Bucharest - Faculty of Applied Chemistry and Material Science,

3 - University of Bucharest - Faculty of Biology,

4 - National Institute for Laser, Plasma and Radiation Physics, Magurele city.


SolariX-XR, QqqFT-ICR HR mass spectrometer with 15 T magnet



EU FT-ICR MS



Interfaces:
MALDI
GC
LC
Direct Probe
Infusion








This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 731077.

Synthesis and characterisation of magnetic core nanoparticles ($\text{Fe}_3\text{O}_4/\text{PABA}$)





EU FT-ICR MS



The synthesis and characterisation of magnetite core-shell nanoparticles with secondary silica shell ($\text{Fe}_3\text{O}_4/\text{PABA}/\text{SiO}_2$)






EU FT-ICR MS



Fe₃O₄-PABA-SiO₂ sp


Absorbance



Dip-coating thin films




EU FT-ICR MS





EU FT-ICR MS


Nanofilms obtained by dip-coating method from dispersions of core-shell magnetite nanoparticles with secondary silica shell ($\text{Fe}_3\text{O}_4\text{-PABA-SiO}_2$) loaded with DPPH radical








EU FT-ICR MS

Nanofilms obtained by dip-coating method from dispersions of core-shell magnetite nanoparticles with secondary silica shell (Fe_3O_4 -PABA- SiO_2) loaded with DPPH radical









EU FT-ICR MS





Synthesis and characterization of MCM-41




MCM - Mesoporous material loaded with Praziquantel



EU FT-ICR MS





a.u. * 1000



Thin films (dip coating) of biopolymer encapsulated MCM41 np's loaded with Harpagophytum extract



EU FT-ICR MS



Higher thickness film



Deposition discontinuity

Low thickness film



ITO/glass slide



MAPLE: MATRIX Assisted Pulsed Laser Evaporation




EU FT-ICR MS






EU FT-ICR MS



a.u. *10e-3




Suport: Si
 Fe_3O_4 / PABA / Ag / PABA / Quercetin
Ion: 269.3892 m/z



Suport: Si
 Fe_3O_4 / PABA / Ag / PABA / Copper phthalocyanine
Ion: 576.6982 m/z



EU FT-ICR MS



Acknowledgment



EU FT-ICR MS

THANK YOU FOR YOUR ATTENTION !

THIS RESEARCH WAS FUNDED BY THE PROJECTS:

- **EUROPEAN REGIONAL DEVELOPMENT FUND THROUGH COMPETITIVENESS OPERATIONAL PROGRAM 2014-2020, PRIORITY AXIS 1, PROJECT NO. P_36_611, MYSMIS CODE 107066, INNOVATIVE TECHNOLOGIES FOR MATERIALS QUALITY ASSURANCE IN HEALTH, ENERGY AND ENVIRONMENTAL - CENTER FOR INNOVATIVE MANUFACTURING SOLUTIONS OF SMART BIOMATERIALS AND BIOMEDICAL SURFACES – INOVABIOMED**
- **„COMPOZIT MULTIFUNCTIONAL PE BAZA DE MATRICE SILICA-ORGANICA TRANSPOZABILA PENTRU INOVATII DE PRODUSE SI FORMULARI PARTICULARIZATE IN INDUSTRIA ALIMENTARA SI FARMACEUTICA”, COD SMIS 122027, CONTRACT NO. 257/09.06.2020.**